Dycal®

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Dycal®

Radiopaque Calcium Hydroxide Composition

INSTRUCTIONS FOR USE – ENGLISH

CAUTION: This is a medical device. For dental use only.
USA: Rx only.

1. PRODUCT DESCRIPTION

Dycal® Calcium Hydroxide Liner is a two-component, rigid-setting, self-curing material designed for use in direct and indirect pulp capping and as a protective liner under dental adhesives, varnishes, filling materials, cements, and other base materials. It will not inhibit the polymerization of acrylic and composite restorations.

1.1 Indications

Dycal® Liner is indicated for:

1. Application to exposed, vital pulp tissue (direct pulp capping).
2. Application to dentin as a protective barrier between restorative materials and deep vital dentin (indirect pulp capping) or where dentin to restorative material contact is not desired.

1.2 Contraindications

None known.

1.3 Delivery forms (some delivery forms may not be available in all countries)

Dycal® Liner is available in:

- collapsible laminate tubes for easy dispensing
- 2 shades: Ivory and dentin

1.4 Composition

- Dycal® Liner Base paste: disalicylate ester of 1,3-butylene glycol; calcium phosphate; calcium tungstate; zinc oxide; iron oxide
- Dycal® Liner Catalyst paste: calcium hydroxide; ethyl toluenesulfonamide; zinc stearate; titanium dioxide; zinc oxide; iron oxide

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapter of these Instructions for Use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

1. Dycal® Liner contains disalicylate ester of 1,3-butylene glycol; contains ethyl toluenesulfonamide; components and mixed material may be irritating to skin, eyes and oral mucosa, and may cause allergic contact dermatitis in susceptible persons.
 - **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
 - **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. Avoid prolonged or repeated exposure to skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
 - **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.
2. Dycal® Liner should not be used with patients who have a history of hypersensitivity to any of the components.

2.2 Precautions

1. This product is intended to be used only as specifically outlined in the Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and sole responsibility of the practitioner.
2. Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
3. The solubility and essentially basic nature of Dycal® Liner composition requires that it be used only in situations where it can be protected adequately from the intraoral environment. After the Dycal® Liner has completely set, cover with cavity varnish, dentin adhesive, or restorative material per manufacturer's instructions. Do not place on enamel, cavity margins or leave exposed to oral environment (see Step-by-Step Instructions for Use).
4. Insufficient data exist to support the use of Dycal® Liner as base, filling material, core build-up or luting cement.
5. Insufficient data exist to support the use of Dycal® Liner as a root canal dressing, medication or repair material.
6. Contact with saliva, blood and sulcus fluid during application may cause failure of the pulp capping procedure. Use adequate isolation such as rubber dam.
7. Pulpal hemorrhage must be controlled prior to the application of Dycal® Liner. If hemorrhage cannot be controlled, consideration should be given for the initiation of endodontic therapy over direct pulp capping procedures.
8. Studies have shown Dycal® Liner to be adequately strong to resist the packing pressures upon restorative materials during placement by usual procedures. Avoid the use of excessively small pluggers and/or high packing pressures, which might cause fractures of the material.
9. Apply Dycal® Liner as a thin layer; use of multiple layers or one thick layer does not improve functionality.
10. Dycal® Liner should be mixed in equal volumes (1:17 to 1:100 by weight) for optimum performance. The mixed material should be homogeneous and streak free prior to application. Variations may affect the radiopacity, strength, and durability of the material.
11. Tooth size, shape, dental radiopacity exposure conditions, and restorative material may affect the ability to radiographically detect the material.
12. Increased humidity and elevated temperature will reduce available work time and set time. Dycal® Liner sets faster in the mouth than on the bench.
13. The Dycal® Liner base and catalyst containers should be tightly closed immediately after use.
14. Do not use instruments which have not been properly (re)processed. Before re-use, instruments must be properly cleaned/decontaminated to reduce the risk of cross contamination.
15. The tubes cannot be reprocessed. To prevent tubes from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that they are handled outside of the patient treatment operator with clean/disinfected gloves.

Interactions

None known.

2.3 Adverse Reactions

Product may irritate the eyes and skin. **Eye contact:** irritation and possible corneal damage.

Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin.

Mucous Membranes: inflammation (see Warnings).

2.4 Storage

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store with original caps tightly closed in a well ventilated place at temperatures between 2 °C and 24 °C/36 °F and 75 °F.
- Allow material to reach room temperature prior to use.
- Protect from moisture.
- Do not freeze.
- Do not use after expiration date.

3. STEP-BY-STEP INSTRUCTIONS

3.1 Direct Pulp Capping

- 3.1.1. Under rubber dam isolation, complete cavity preparation outline using high-speed burs under constant water-cooling.
- 3.1.2. If caries is present, completely excavate using low-speed and/or hand instrumentation.
- 3.1.3. Rinse the cavity and exposure site(s) with 2.6%-5% NaOCl. Heavy bleeding may be controlled with a cotton pellet moistened with sterile saline.
- 3.1.4. Gently dry preparation with cotton pellet. Avoid desiccation.
- 3.1.5. Dispense equal volumes of base and catalyst pastes on the parchment paper pad provided. Replace container caps. Using a Dycal® Liner applicator, stir immediately to mix thoroughly until a uniform color is achieved. Do not over-spatulate. Complete mixing within 10 seconds. **NOTE:** Do not attempt to control the setting time by increasing or decreasing the amount of catalyst dispensed.
- 3.1.6. Using the ball-pointed Dycal® Liner applicator or similar instrument, place the mix directly on the exposed pulp and cavity dentin judged to be less than 1.0 mm remaining thickness in a thin layer. Avoid placing Dycal® Liner on enamel or the margins of the cavity. Avoid placing a large bulk of material. Material thickness should be approximately 0.8 mm-1 mm.
- 3.1.7. Allow the Dycal® Liner to completely set. The mixed material will set in approximately 2-3 minutes on the mixing pad under normal room conditions (21 °C/70 °F with 50% relative humidity). Set time is shorter in the mouth due to moisture and temperature.
- 3.1.8. Remove any set excess from retention areas, enamel, and/or margins with a sharp spoon excavator or a bur.
- 3.1.9. Place desired adhesive, base, and/or restoration following manufacturer's instructions.
- 3.1.10. At the next appointment, assess the pulp vitality. Pulp vitality and status should be assessed radiographically every three to six months or as needed.

3.2 Indirect pulp capping, protective barrier

- 3.2.1. Under suitable isolation, complete cavity preparation and caries removal. Wash cavity thoroughly with water spray and air dry.
- 3.2.2. Dispense and mix Dycal® Liner components as outlined above in section 3.1.5.
- 3.2.3. Apply mixed material to desired dentin surfaces as outlined above in sections 3.1.6.-3.1.8. Quick application of the mixed material into the cavity takes advantage of the essential flowability of the mass for efficient placement before setting starts. Completely covering all cavity dentin as desired may require carrying additional material from pad to dentin surfaces. **NOTE:** If subsequent use of a dentin bonding agent is desired, place Dycal® Liner only on the deepest (less than 1 mm remaining) dentin, leaving the rest of the cavity surface free for bonding.
- 3.2.4. Following completion of set (see section 3.1.7.) remove any set excess from retention areas, enamel, and/or margins with a sharp spoon excavator or a bur.
- 3.2.5. Complete restoration per restorative material manufacturer's instructions for use.

4. HYGIENE, PROCESSING AND DISPOSAL

Cross-contamination – To reduce risk of Infection

- Reprocess reusable products as described below.
- Tubes cannot be reprocessed. Dispose of contaminated tubes in accordance with local regulations.
- When not in use, all supplies and instruments should be in covered storage, such as drawers and cabinets, and away from potential contamination.
- During treatment clinicians with patient contact should not handle the tubes.
- To prevent the tubes from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the tubes are handled with clean/disinfected gloves, outside of the patient treatment operator.
- Dispensing materials with clean/disinfected gloves or hands in a separate room, bringing only what will be used into the operator is strongly recommended.
- Do not reuse tubes if contaminated.

To prevent tubes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended.

The use of protective barriers is an additional precautionary measure against gross debris but not against all contamination. The barrier is designed for single use and must be disposed of after each use in accordance with local regulations.

4.1 Dycal® Placement Instrument and Tubes

Process		✓ Acceptable
Cleaning	Mechanical (Ultrasonic)	✓
	Manual	✓
	Automated (Washer/Disinfector)	✓
Disinfection	Automated (Washer/Disinfector)	✓
	Manual	✓
	Immersion	
Sterilization	Steam Autoclave	✓
	Immersion	
Single Use	Dispose of after Use	
Reusable/Not Reprocessable (Tubes)	Dispose of if contaminated	✓

Instructions for cleaning, disinfecting and sterilizing Dycal® Placement Instrument	
Warnings	- High-level disinfection has not been validated as a terminal process for the Dycal® placement instrument. Disinfection is not recommended as a terminal process. Proceed to Sterilization following cleaning and any optional/intermediate disinfection process. - Steam autoclaving sterilization is appropriate and recommended for the Dycal® placement instrument. - Follow proper infection prevention activities, such as proper hand washing and donning new puncture and chemical gloves at pertinent steps. - Do not use hard wire brushes to clean the Dycal® placement instrument as scratching and mechanical degradation can occur. - Use only a disinfecting solution which is approved for its efficacy, EPA registered (and/or Health Canada approved) and use in accordance with the IFU of the disinfecting solution manufacturer. - Do not allow the Dycal® placement instrument to exceed 137 °C/279 °F. - Do not use phenol-based glutaraldehyde. - Always use a pH neutral instrument cleaning solution if authorized. - Sterilize metal and plastic components in separate pouches to avoid damage to components. - If the integrity of the pouch has been compromised, the device should be processed again before use. - All devices should be cleaned and sterilized before each use. For the countries that require a three-step process, in which disinfection is required prior to sterilization, all devices must be cleaned, disinfected, and sterilized before use.
Limitations on Reprocessing	- Repeated reprocessing has minimum effect on these instruments. End of life is normally determined by wear and damage due to use. - Cold liquid immersion disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.
Initial Treatment at the Point of Use	- Do not allow residue or any form of contamination to dry on the device. Remove excess soil/gross debris by wiping, brushing, and/or rinsing under water. - Transport to processing area and follow validated processing instructions below. - It is recommended that the device be processed as soon as is reasonably practical following use.
Preparation before cleaning	No particular requirements.
Cleaning: Mechanical	- The chemical manufacturer's instructions for the cleaning and disinfection solution should be followed, observing concentration rates, contact times, fill level, degassing processes, etc. - The recommended solution is Resurge® Instrument Cleaning Solution mixed to 0.5 oz/gal (14.8 ml/3.8 l), at 35 °C-40 °C/95 °F-104 °F, FD 370 Cleaner, Dürr Dental AG, Germany or equivalent enzymatic cleaning solution. - Immerse Dycal® placement instrument in an ultrasonic bath containing the prepared solution for at least 20 minutes unless otherwise specified by the time specified by the solution manufacturer. - Rinse under tap water of at least drinking quality for a minimum of 30 seconds. - Allow to air dry. - Visually inspect for visible soil, and repeat the cleaning and disinfection procedure if necessary.
Cleaning: Manual	- Clean Dycal® placement instrument with a soft brush, tap water, and liquid soap until free of visible contamination. - Remove cleaning solution residue with a cloth soaked with tap water or rinse under warm tap water for 2 minutes. - Dry Dycal® placement instrument with a lint-free, single-use cloth or allow to air dry in a clean area. - Visually inspect for visible soil, and repeat the cleaning procedure if necessary.
Disinfection: Manual (Intermediate only, not terminal process)	Disinfection by wiping. - No manual or automated terminal disinfection process has been validated. Device has been shown to be compatible and disinfected with an alcohol-based, tuberculocidal, quaternary ammonium solution, approved according to local regulations, and used according to disinfectant solution manufacturer's Instruction for Use. - After the desired contact time, remove disinfectant residue with a single-use cloth soaked in deionized water. - Dry the Dycal® placement instrument using a dry, lint-free disposable cloth or disposable wipe. - Proceed to Sterilization following manual cleaning and any additional disinfection process.
Cleaning and disinfection: Automated	- A thermal disinfectant that complies with ISO 15883 may be used for automatic cleaning and disinfection. - Use only properly maintained, calibrated, and approved washer disinfectors according to ISO 15883. - Follow the manufacturer's instructions to select the solutions and programs. - Remove excess soil with a disposable cloth/paper wipe. - In case of gross contamination of Dycal® placement instrument, remove it using a soft brush and running water. - Place Dycal® placement instrument in the washer-disinfector allowing water and detergent to enter into and drain out of device openings. - Run washer-disinfection program with A0 value ≥ 3000 (e.g. 5 min at ≥ 90 °C/194 °F) (e.g. Miele Vario TD) using appropriate detergents (e.g. neodisher® MediClean and neodisher® Z, both from Dr. Weigert, Germany or equivalent). - After the cycle is complete, visually inspect the Dycal® placement instrument for degradation due to extended use and reprocessing. If the device is discolored, cracked, worn, distorted, etc., discard and do not use. - Proceed to Sterilization following Automated Cleaning and Disinfection.
Drying	The drying instructions are incorporated in the cleaning and disinfection section above.

Maintenance, Inspection and Testing	- Visually inspect the Dycal® placement instrument for degradation due to extended use and reprocessing. If the device is discolored, cracked, worn, distorted, etc., discard and do not use. - Dycal® placement instrument requires no additional maintenance and should not be oiled.
Packaging	- Ensure that the Dycal® placement instrument is completely dry before packaging for sterilization. - A lint-free, disposable cloth may be used to dry Dycal® placement instrument. - Place each Dycal® placement instrument in a separate appropriately sized steam-sterilization pouch. - Uses of FDA-cleared, ISO 11607 compliant Paper/plastic steam sterilization pouches (e.g., AssurePlus® Sterilization Pouches) are recommended.
Sterilization/ Steam autoclaving	- The Dycal® placement instrument must be sterilized utilizing steam autoclave sterilization. NOTE: Dycal® placement instruments may be sterilized by Dynamic-Air-Removal/Pre-Vacuum steam autoclave. The following Pre-Vacuum cycles may be used: 134 °C/273 °F for 3 minutes 30 seconds (recommended only for non-US markets). NOTE: Dycal® placement instruments may be sterilized by Gravity Displacement steam autoclave. The following Gravity cycles may be used: 135 °C/275 °F for 10 minutes with a minimum drying time of 30 minutes. An alternate method of sterilization is to place non-bagged instruments into the steam sterilizer and run at one of the above listed cycles.
Storage	- Do not remove from the pouch until ready for use to prevent contamination. - Instruments sterilized non-bagged should be used immediately. - Store at room temperature, away from moisture or excessive humidity. - To prevent contamination, store processed device in covered storage such as a drawer or cabinet until use. - Prior to use, inspect the pouch. If the integrity of the pouch has been compromised, the device should be processed again before use.
Additional Information	- Prior to use, inspect the device. Discard any Dycal® placement instrument which has become damaged, worn, or distorted. - The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for use/reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and trained personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences. Use of other reprocessing procedures/methods is at the discretion and sole responsibility of the practitioner.
Manufacturer Contact	Within the United States, call Dentsply Sirona at 1-302-422-4511. For areas outside the United States, contact your local Dentsply Sirona representative.

4.2 Disposal

- Dispose of in accordance with local regulations.
- Dispose of disinfected mixing spatula, tubes, and packages in normal waste stream in accordance with local regulations.
- See SDS for further details.

5. LOT NUMBER, EXPIRATION DATE AND CORRESPONDENCE

- Do not use after expiration date. ISO standard uses: "YYYY-MM-DD."
- The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number on package
 - Expiration date
- Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
- For the summary of safety and clinical performance for this product please visit <https://ec.europa.eu/tools/eudamed> and/or <https://www.dentsplysirona.com/en/service-contact/download-center.html>. The Device Identification: Basic UDI-DI ++D002BASE55 can be used to search.

163-DVO

© 2025 Dentsply Sirona Inc. All Rights Reserved. **523001WEB** (R 1/9/25)

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